



mim


Enterprise

Solutions

# Why MIM?

MIM provides integral solutions that organize and improve department workflows for centers with diverse software and equipment. Through constant innovation, MIM continues to develop technologies that exceed the expectations of large institutions.

- **Flexibility** - By supporting virtually all OEM vendors and third-party systems, MIM is the solution that fits in every environment
- **Reading** - Review automatically linked images from multiple vendors simultaneously
- **Fusion** - Fuse any combination of images from any vendor including CT, MR, PET, SPECT, US, XA and most 2D images
- **Storage** - Expandable storage to accommodate both DICOM and non-DICOM objects such as PDF, JPG, XLS, and more
- **Integration** - MIM supports integrations with many major vendors, and our API framework is simple to implement for your RIS, EMR, PACS, etc., to further optimize workflow and accessibility




**Vendor Neutral:** MIM integrates with virtually all OEM vendors and third-party systems.

## mim Zero Footprint™

Access MIM Anywhere and Everywhere.

- No local application servers
- Access on any device with a web browser
- Delivered via encrypted HTML5 virtualization
- Automated mirroring - local database and cloud database are synced
- Allows for transparent user experience regardless of access point



# mim Assistant™

The Engine of Automation.

## Clinical Workflow Management

- Automatically process incoming patient data through a rule-based environment
- Automatically run MIM Workflows™
- Auto-contour CT's and save the structures to the patient list
- Automatically fuse images before you arrive at your workstation to review
- Auto-generate plan reports in Plan Review™

## Data Management

- Auto-sync and compress data to MIMcloud™ for quick access
- Auto-query diagnostic images
- Auto-route studies and automate fetching protocols

## User Management

- Customize preferences, workflows, and privileges for all users
- LDAP/AD integration allows for the use of existing permissions, profiles, and login credentials
- Task-list functionality allows for assigning studies or statuses to coordinate the planning process

## Batch Processing

- Large populations of patient data can be processed using custom workflows
- Most functions in MIM can be automated, so extracting information, organizing, and feeding reports or spreadsheets is quick, reducing errors on large projects

# mim Workflows

The Framework of Automation.

- Automate complex or time-consuming tasks like contouring, TV expansions, and deformable image registration
- Create flexible hanging protocols for more efficient reading including serial study review, HN & WB, and WB with Legs with automated linking
- Customization for full control of automated processes
- Drag-and-drop workflow design interface – no coding involved

## Session Save

Preserve Your Work and Instantly Resume Where You Left Off.

- Organizes and expedites patient management by allowing clinicians to preserve complex processes for later review
- Easily access patient data for tumor board/chart rounds
- Hand off a nuclear medicine case from a technologist to a physician, exactly as it was prepared



A Clinical System that  
Excels at Research.

#### Native support for MATLAB® and Java®

- Write your own code and utilize it within MIM
- Leverage MIM Extensions™ and MIM Assistant for batch processing

#### Expedite Research

- MIM's intuitive interface lets you focus on your specific project and avoid redundant tasks such as building a UI, library, or tools

#### Translational Research

- Immediately apply lab developments in the clinic

## User Preference Management and Security

- User access auditing
- Permission based access to PHI
- Easy management for administrators for both preferences and security settings

## Support

### A Customer-centric Approach.

- Clinical Support Engineers directly available via the web or by phone for immediate consultation
- Dedicated Support Manager assigned to your site responsible for ongoing training and routine follow-up
- Customer feedback is highly-valued and documented to continually improve our product
- Pager system for quick call-backs on nights and weekends

To learn more, call **866-421-2536** or visit **[mimsoftware.com/contact](http://mimsoftware.com/contact)** to schedule a presentation of MIM.



#### **MIM SOFTWARE INC.**

25800 Science Park Drive - Suite 180  
Cleveland, OH 44122

866-421-2536

[www.mimsoftware.com](http://www.mimsoftware.com)

[info@mimsoftware.com](mailto:info@mimsoftware.com)

#### **MIM SOFTWARE**

Beijing Co., Ltd.

北京明维视景医疗软件开发有限公司  
地址: 北京市海淀区学院路51号首享科  
技大厦809室  
邮编 100191

电话 86-10-82626960

邮箱 [info@mimsoftware.com](mailto:info@mimsoftware.com)



#### **MIM SOFTWARE**

Brussels BV

Drukpersstraat 4  
1000 Brussel  
Belgium

© MIM Software Inc.

2024 All Rights Reserved

**TD-442**

18 Nov 2024



#### **Emergo Europe**

Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands



#### **MedEnvoy Switzerland**

Gotthardstrasse 28  
6302 Zug  
Switzerland

**CE**  
**1639**